

# How to Update my Guide

- 1. Logon to the server.
- 2. Open setupmgr.exe (Note This location will depend on the site).

#### 3. This screen will open.

			Sector Sector				
Instance	Product Name	Component Name	Install/Repair T	his Vers	iion	Config	Uninst
PROD APP	Cardiac Care	Cardiac Care App Server	10.0.2044.10006	•	0	*	×
PROD INTERFACE	Cardiac Care	Interface Engine	10.0.2044.10006	•	0	*	×
PROD LOGIN	Cardiac Care	Windows Domain Login	10.0.2044.10006	•	0	*	×
PROD WEB	Cardiac Care	Cardiac Care	10.0.2044.10006	•	0	*	×
TEST APP	Cardiac Care	Cardiac Care App Server	10.0.2044.10006	•	0	*	X
TEST INTERFACE	Cardiac Care	Interface Engine	10.0.2044.10006	•	C	*	×
TEST LOGIN	Cardiac Care	Windows Domain Login	10.0.2044.10006	•	0	*	×
TEST WEB	Cardiac Care	Cardiac Care	10.0.2044.10006	•	Ð	*	×
continue 💌	Bulk Upgrade	Create New Instance		Remove Unused Ver		sions	
Update Cardiac Care to 10.0.2044.10007				Stage		1	
-							
When install	ing/repairing instances, allow datab	base operations to take place				CL	

### 4. Click on the Update Bar.

Update Cardiac Care to 10.0.2044.10007

#### 5. A progress bar will appear.

Downloading Update



## 6. Select All.

5.

🍥 Cedaron Form Builder Setup		×
PROD APP PROD INTERFACE PROD LOGIN PROD WEB TEST INTERFACE TEST INTERFACE TEST LOGIN TEST WEB	10.02044.10006 10.02044.10006 10.02044.10006 10.02044.10006 10.02044.10006 10.02044.10006 10.02044.10006 10.02044.10006	Select All Select None Do Sites Sequentially
< Prev ☑When in ☑When in	stalling/repairing instances, allow database operations to take place stalling/repairing instances, allow database downgrade operations to take place	Next->
Select Next		