


Updating Cento Flash Reports

Cedaron periodically adds and updates the built-in Cento flash reports. When new flash reports are released by Cedaron, a CardiacCare admin will have to log into Cento Analytics and update the flash reports in the organization's Cento Analytics instance.




How to Download/Update Flash Reports

Follow these steps to download/update flash reports:


1. Launch **Cento Analytics**. This is done in CardiacCare by hovering over **Utilities** in the menu bar, then selecting **Analytics**.




2. Go to **Administration**. Can click on either Administration, in the Cento menu bar or the Cento home screen dashboard (both highlighted in red below).



3. Go to **Manage Flash Reports** (shown in red below).



4. Check the boxes next to each registry name then click on **Click Here to Download Flash Reports**.



- Once the flash reports are updated, the All Report Templates page will load and there will be a banner notification stating **Flash Reports Updated Successfully**.

The screenshot shows the Cedaron web application interface. At the top, a navigation bar includes 'Cenlo', 'All Report Templates', 'New Report', and 'Administration'. On the right, it shows the user 'TrainingAuthVM/cm_ebckr' and a 'Log Off' link. A blue banner notification at the top center reads 'Flash Reports Updated Successfully'. Below this, the page title is 'All Report Templates'. There are three tabs: 'My Reports', 'Shared Custom Reports', and 'Flash Reports'. The 'Flash Reports' tab is active. The main content area features a table of report templates with columns for 'Report Template', 'Report Type', 'Last Run', 'Tags', and 'Run Report'. The table lists 10 report templates, including STS 2.9 Clinical Dashboards and CathPCI v5 reports. A pagination bar at the bottom indicates 'Showing 1 to 10 of 208 records' and includes page numbers 1 through 21.

| Report Template | Report Type | Last Run | Tags | Run Report |
|---|---------------|--------------------|-----------------------|------------|
| STS 2.9 Clinical Dashboard - Isolated CABG - Process of Care | Summary | 4/19/2021 4:40 PM | STS Adult Cardiac 2.9 | Run Report |
| STS 2.9 Clinical Dashboard - Isolated CABG - Outcomes - 1 | Summary | 4/28/2021 4:03 PM | STS Adult Cardiac 2.9 | Run Report |
| STS 2.9 Clinical Dashboard - Isolated CABG - Outcomes - 2 | Summary | 4/19/2021 7:27 PM | STS Adult Cardiac 2.9 | Run Report |
| STS 2.9 Clinical Dashboard - Isolated AVR - Process of Care | Summary | 10/7/2020 10:41 PM | STS Adult Cardiac 2.9 | Run Report |
| STS 2.9 Clinical Dashboard - Isolated AVR - Outcomes - 1 | Summary | 10/7/2020 10:41 PM | STS Adult Cardiac 2.9 | Run Report |
| STS 2.9 Clinical Dashboard - Isolated AVR - Outcomes - 2 | Summary | 10/7/2020 10:42 PM | STS Adult Cardiac 2.9 | Run Report |
| CathPCI v5 - Case Volumes | Pivot Summary | 5/6/2021 5:22 PM | ACC/NCDR CathPCI v5.0 | Run Report |
| CathPCI v5 - Metric 2: PCI Procedures with Positive Stress or Imaging Study | Pivot Summary | 5/6/2021 6:48 PM | ACC/NCDR CathPCI v5.0 | Run Report |
| CathPCI v5 - Metric 9: P2Y12 inhibitor prescribed at discharge | Pivot Summary | 5/6/2021 6:22 PM | ACC/NCDR CathPCI v5.0 | Run Report |
| CathPCI v5 - Metric 10: Statin prescribed at discharge | Pivot Summary | 3/4/2021 6:13 PM | ACC/NCDR CathPCI v5.0 | Run Report |